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NUROTRON DIGEST

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文献摘要.....	1
Stable benefits of bilateral over unilateral cochlear implantation after two years: A randomized controlled trial.	1
术后两年双侧优于单侧人工耳蜗植入——双侧稳定受益：随机对照试验 ...	2
Bilateral cochlear implantation in children: a systematic review and best-evidence synthesis.	3
儿童双侧人工耳蜗植入：系统综述与最佳证据整合	4
Sequential bilateral cochlear implantation in the adolescent population.....	5
青少年人群中相继双侧人工耳蜗植入术	6
Bilateral cochlear implantation for hearing-impaired children: criterion of candidacy derived from an observational study.	7
听障儿童双侧人工耳蜗植入：基于观察研究的候选标准	8
Bilateral cochlear implantation.	9
双侧人工耳蜗植入术	10
The effect of device use after sequential bilateral cochlear implantation in children: An electrophysiological approach.	11
儿童双侧相继人工耳蜗植入的术后设备使用效果：电生理学方法	12
Bilateral sequential adult cochlear implantation: who should receive priority in the context of a constrained health care system	13
成人双侧相继人工耳蜗植入：在有限的卫生保健系统中谁应该获得优先救助？	14
人工耳蜗双侧植入的优势.....	15
人工耳蜗植入的双耳聆听研究进展	15
双侧人工耳蜗植入者在噪声环境下的言语辨别能力.....	16
双侧同期人工耳蜗植入效果分析	16
双侧人工耳蜗植入术后言语识别效果的评估	17
双侧人工耳蜗植入对儿童听觉识别能力影响的 个案研究.....	17

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Stable benefits of bilateral over unilateral cochlear implantation after two years: A randomized controlled trial.

Van Zon A¹, Smulders YE¹, Stegeman I¹, Ramakers GG¹, Kraaijenga VJ¹, Koenraads SP¹, Zanten GA¹, Rinia AB², Stokroos RJ³, Free RH⁴, Frijns JH⁵, Huinck WJ⁶, Mylanus EA⁶, Tange RA¹, Smit AL¹, Thomeer HG¹, Topsakal V¹, Grolman W¹

1. Brain Center Rudolf Magnus, University Medical Center Utrecht, Utrecht, the Netherlands.
2. Department of Otorhinolaryngology-Head and Neck Surgery, University Medical Center Utrecht, Utrecht, the Netherlands.
3. Department of Otorhinolaryngology-Head and Neck Surgery, Maastricht University Medical Center, Maastricht, the Netherlands.
4. Graduate School of Medical Sciences, Research School of Behavioural and Cognitive Neurosciences, University Medical Center Groningen, Groningen, the Netherlands.
5. Leiden Institute for Brain and Cognition, Leiden University Medical Center, Leiden, the Netherlands.
6. Donders Institute for Brain, Cognition and Behaviour, Radboud University Medical Center, Nijmegen, the Netherlands.

Abstract

To investigate hearing capabilities and self-reported benefits of simultaneous bilateral cochlear implantation (BiCI) compared with unilateral cochlear implantation (UCI) after a 2-year follow-up and to evaluate the learning effect of cochlear implantees over time. Multicenter randomized controlled trial. Thirty-eight postlingually deafened adults were included in this study and randomly allocated to either UCI or simultaneous BiCI. Our primary outcome was speech intelligibility in noise, with speech and noise coming from straight ahead (Utrecht-Sentence Test with Adaptive Randomized Roving levels). Secondary outcomes were speech intelligibility in noise with spatially separated sources, speech intelligibility in silence (Dutch phoneme test), localization capabilities and self-reported benefits assessed with different quality of hearing and quality of life (QoL) questionnaires. This article describes the results after 2 years of follow-up. We found comparable results for the UCI and simultaneous BiCI group, when speech and noise were both presented from straight ahead. Patients in the BiCI group performed significantly better than patients in the UCI group, when speech and noise came from different directions ($P = 0.01$). Furthermore, their localization capabilities were significantly better. These results were consistent with patients' self-reported hearing capabilities, but not with the questionnaires regarding QoL. We found no significant differences on any of the subjective and objective reported outcomes between the 1-year and 2-year follow-up. This study demonstrates important benefits of simultaneous BiCI compared with UCI that remain stable over time. Bilaterally implanted patients benefit significantly in difficult everyday listening situations such as when speech and noise come from different directions. Furthermore, bilaterally implanted patients are able to localize sounds, which is impossible for unilaterally implanted patients. .

Article source:

Van Z A, Smulders Y E, Stegeman I, et al. Stable benefits of bilateral over unilateral cochlear implantation after two years: A randomized controlled trial.[J]. Laryngoscope, 2017, 127(5):1161.

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Van Z A, Smulders Y E, Stegeman I, et al. Stable benefits of bilateral over unilateral cochlear implantation after two years: A randomized controlled trial.[J]. Laryngoscope, 2017, 127(5):1161.

Bilateral cochlear implantation in children: a systematic review and best-evidence synthesis.

Lammers MJ¹, van der Heijden GJ, Pourier VE, Grolman W.

Department of Otorhinolaryngology-Head and Neck Surgery, University Medical Center Utrecht, Utrecht; Brain Center Rudolf Magnus, University Medical Center Utrecht, Utrecht.

Abstract

Objectives/hypothesis: To evaluate the effectiveness of bilateral cochlear implantation over unilateral implantation in children with sensorineural hearing loss.

Data sources: Pubmed, Embase, and Web of Science.

Review methods: All studies comparing a bilateral cochlear implant group with a unilateral implant group were included.

Results: Twenty-one studies compared a bilateral cochlear implant group with a unilateral group. No randomized trials were identified. Due to the clinical heterogeneity, statistical pooling was not feasible and a best-evidence synthesis was performed. The results of this best-evidence synthesis indicate the positive effect of the second implant for especially sound localization and possibly for preverbal communication and language development. There was insufficient evidence to make a valid comparison between bilateral implantation and a bimodal fitting.

Conclusion: Although randomized trials are lacking, the results of our best-evidence synthesis indicate that the second cochlear implant might be especially useful in sound localization and possibly also in language development.

Key Words: Cochlear implantation; best-evidence synthesis; bilateral; bimodal; deafness; hearing loss; systematic review; unilateral implantation; best-evidence synthesis; bilateral; bimodal; deafness; hearing loss; systematic review; unilateral.

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Sequential bilateral cochlear implantation in the adolescent population.

Friedmann DR¹, Green J¹, Fang Y², Ensor K¹, Roland JT¹, Waltzman SB¹.

1. Department of Otolaryngology-Head & Neck Surgery, New York University School of Medicine, New York, New York, U.S.A.

2. Division of Biostatistics, Department of Population Health, New York University School of Medicine, New York, New York, U.S.A.

Abstract

Objectives: To examine the variables affecting outcomes for sequential bilateral cochlear implantation patients in the adolescent population.

Study design: Retrospective chart review at tertiary care center.

Methods: Main outcome measures were open set speech recognition tests at the word (Consonant-Nucleus-Consonant/Phonetically Balanced Kindergarten List Test [CNC/PBK]) and sentence levels in noise (Hearing in Noise Test-Noise [HINT-N]) in different test conditions with respect to the age at first and sequential implantation, as well as the interval between implants.

Results: Despite a mean age at sequential implantation of 13.5 years, sequential bilateral implanted adolescents revealed significant improvement in the sequential cochlear implant (CI2) ear. The mean time interval between implants was 8.2 years. A wide range of performance was noted, and age at implantation and interval between first cochlear implant (CI1) and CI2 did not predict outcome. Mean CNC/PBK score with CI1 alone was 83.0%, with the CI2 alone was 56.5%, and with bilateral implants was 86.8%. Sentence scores (HINT-N) were 89.5% for CI1, 74.2% for CI2, and 94.4% for bilateral CI condition. The clinical relevance of these enhanced perception abilities requires attention to individual device use, performance with the first implant, and subjective benefits reported by patients.

Conclusions: Bilateral sequential cochlear implantation leads to improved speech perception in the adolescent population and should be considered in this population, even after a long period of deafness and despite a prolonged interval between implants. Numerous factors affect the ability to predict performance, but age at implantation and interimplant interval were not correlated with performance measures. Extensive preoperative counseling and individualized evaluation are critical to ensure that patients and families understand the range of possible outcomes.

Key Words: Adolescent, pediatric cochlear implants; bilateral; sequential.

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Bilateral cochlear implantation for hearing-impaired children: criterion of candidacy derived from an observational study.

Lovett RE¹, Vickers DA, Summerfield AQ.

¹Ear Institute, University College London, London, United Kingdom; and ²Department of Psychology and Hull York Medical School, University of York, York, United Kingdom.

Abstract

Objectives: Policy-makers have struggled to define the minimum degree of hearing impairment at which children should be offered cochlear implants rather than the less invasive alternative of acoustic hearing aids. This study compared outcomes for children with bilateral cochlear implants and children with bilateral hearing aids, to determine a criterion of candidacy for pediatric bilateral cochlear implantation.

Design: This observational study measured the listening skills of children who received routine audiological care in the United Kingdom. Participants were recruited from hospitals, educational services, and charities. Eligibility criteria included a diagnosis of hearing impairment before 31 months of age and pure-tone thresholds greater than or equal to 50 dB HL at 2 and 4 kHz bilaterally. Seventy-one children participated, aged 46 to 86 months (mean 64 months). Twenty-eight children used bilateral implants provided in a simultaneous surgery; 43 used bilateral digital hearing aids. The two groups of children were demographically similar in variables that predict outcomes for children with hearing impairment. Children's ability to understand speech was measured using closed-set tests of word discrimination in three conditions: in quiet, in pink noise, and in two-talker babble. For each listening test, an actuarial method was used to compare the distribution of scores from children with cochlear implants and children with hearing aids. The aim was to calculate the unaided pure-tone average (PTA) hearing level at which a child has odds of 4:1 of a better outcome with implants than with hearing aids. The PTA associated with odds of 4:1 has been used previously to define criteria of candidacy for implantation. The main analyses used a four-frequency PTA (mean of unaided thresholds at 0.5, 1, 2, and 4 kHz in the better-hearing ear). Additional analyses used a three-frequency PTA (0.5, 1, and 2 kHz) and two-frequency PTA (2 and 4 kHz).

Results: Odds of 4:1 of a better outcome with implants were associated with a four-frequency PTA of 79, 86, and 76 dB HL for tests of word discrimination in quiet, noise, and babble, respectively. The mean of these three estimates is 80 dB HL. It can be difficult to measure a four-frequency PTA in young children, but a two-frequency PTA typically can be measured. Odds of 4:1 were associated with a two-frequency PTA of 83, 92, and 80 dB HL for tests of word discrimination in quiet, noise, and babble, respectively. The mean of these three estimates is 85 dB HL.

Conclusions: Children with an unaided four-frequency PTA of 80 dB HL or poorer in both ears should be considered candidates for bilateral cochlear implantation. In cases where a four-frequency PTA cannot be measured, the criterion of candidacy should be a two-frequency PTA of 85 dB HL or poorer in both ears. If adopted by policy-makers, these recommendations would expand the provision of cochlear implants among children in England and Wales.

Article source:

Lovett R E, Vickers D A, Summerfield A Q. Bilateral cochlear implantation for hearing-impaired children: criterion of candidacy derived from an observational study.[J]. Ear & Hearing, 2015, 36(1):14-23.

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Lovett R E, Vickers D A, Summerfield A Q. Bilateral cochlear implantation for hearing-impaired children: criterion of candidacy derived from an observational study.[J]. Ear & Hearing, 2015, 36(1):14-23.

Bilateral cochlear implantation.

Wanna GB¹, Gifford RH, McRackan TR, Rivas A, Haynes DS.

Otolaryngologic Clinics of North America.

Abstract

Cochlear implantation (CI) is the standard of care for the treatment of children and adults with bilateral severe-to-profound sensorineural hearing loss. Because the ultimate and continuous goal of CI teams is to improve patient performance, a potential method is bilateral CI. The potential benefits of bilateral CI include binaural summation, squelch, equivalent head shadow for each ear, improved hearing in noise, sound localization ability, and spatial release from masking. The potential disadvantages include additional or prolonged surgical procedure, unproven cost/benefit profile, and the elimination of the ability to use future technologies and/or medical therapies in the implanted ear.

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The effect of device use after sequential bilateral cochlear implantation in children: An electrophysiological approach.

Sparreboom M¹, Beynon AJ², Snik AF², Mylanus EA².

¹ Donders Institute for Brain, Cognition and Behaviour, Department of Otorhinolaryngology, Head and Neck Surgery, Hearing and Implants, Radboud University Medical Centre Nijmegen, The Netherlands. Electronic address: Marloes.Sparreboom@Radboudumc.nl.

² Donders Institute for Brain, Cognition and Behaviour, Department of Otorhinolaryngology, Head and Neck Surgery, Hearing and Implants, Radboud University Medical Centre Nijmegen, The Netherlands.

Abstract

Objectives: In many studies evaluating the effect of sequential bilateral cochlear implantation in congenitally deaf children, device use is not taken into account. In this study, however, device use was analyzed in relation to auditory brainstem maturation and speech recognition, which were measured in children with early-onset deafness, 5-6 years after bilateral cochlear implantation. We hypothesized that auditory brainstem maturation is mostly functionally driven by auditory stimulation and is therefore influenced by device use and not mainly by inter-implant delay.

Methods: Twenty-one children participated and had inter-implant delays between 1.2 and 7.2 years. The electrically-evoked auditory brainstem response was measured for both implants separately. The difference in interaural wave V latency and speech recognition between both implants were used in the analyses. Device use was measured with a Likert scale.

Results: Results showed that the less the second device is used, the larger the difference in interaural wave V latencies is, which consequently leads to larger differences in interaural speech recognition.

Conclusions: In children with early-onset deafness, after various periods of unilateral deprivation, full-time device use can lead to similar auditory brainstem responses and speech recognition between both ears. Therefore, device use should be considered as a relevant factor contributing to outcomes after sequential bilateral cochlear implantation. These results are indicative for a longer window between implantations in children with early-onset deafness to obtain symmetrical auditory pathway maturation than is mentioned in the literature. Results, however, must be interpreted as preliminary findings as actual device use with data logging was not yet available at the time of the study.

Key Words: Bilateral cochlear implants; Device use; Electrically evoked auditory brainstem response; Pediatric; Speech recognition

Article source:

Sparreboom M, Beynon A J, Snik A F, et al. The effect of device use after sequential bilateral cochlear implantation in children: An electrophysiological approach[J]. International Journal of Pediatric Otorhinolaryngology, 2016, 86:161-166.

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Sparreboom M, Beynon A J, Snik A F, et al. The effect of device use after sequential bilateral cochlear implantation in children: An electrophysiological approach[J]. International Journal of Pediatric Otorhinolaryngology, 2016, 86:161-166.

Bilateral sequential adult cochlear implantation: who should receive priority in the context of a constrained health care system?

Forzley B¹, Chen J, Nedzelski J, Lin V, Shipp D, Godlovitch G, Hebert P, Hochman J.

1. Penticton Regional Hospital Renal Program Clinical Instructor, Department of Medicine, University of British Columbia Specialist Champion, British Columbia, Penticton, British Columbia.

Abstract

Resource allocation decisions have become increasingly necessary as the cost of health care habitually increases. Bilateral (second side) adult cochlear implantation (CI) is an example of a novel technology with accruing evidence of benefit, yet expense has limited universal employ. Currently at our centers, bilateral implantation is only provided under research protocol. In this article, we discuss the need for a principled approach concerning the distribution of a second device, both during this period of investigation and if ultimately an insured service. Allocation strategies, while extensively addressed in some arenas, have yet to be developed for second-side sequential adult CI. We advocate that physicians must assume an explicit role when both caring for individual patients as well as administering health care programs. We review social justice theories that inform resource allocation macrodecisions, and include a defence of age-based considerations. Our approach to patient selection for adult second-side CI sequentially considers clinical criteria (directly addressed in the article), a willingness to participate in rigorous research, and a 65 year cut-off. Ultimately, we employ random blinded selection for allocating bilateral CI among the remaining similarly situated individuals. This approach functions impartially and in a manner that is transparent for both patient and physician.

Article source:

Forzley B, Chen J, Nedzelski J, et al. Bilateral sequential adult cochlear implantation: Who should receive priority in the context of a constrained health care system?[J]. *Laryngoscope*, 2013, 123(12):3137-40.

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Forzley B, Chen J, Nedzelski J, et al. Bilateral sequential adult cochlear implantation: Who should receive priority in the context of a constrained health care system?[J]. Laryngoscope, 2013, 123(12):3137-40.

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